

# DISPOSABLE RESPIRATORY DEVICES

## Catalog



**PEPPER MEDICAL**

1-800-647-0172 | [PepperMedical.com](http://PepperMedical.com)

# Contents

Page 2 | **Respiratory Masks**

Page 3 | **Nebulizer Kits**

Page 4 | **Oxygen Tubing**

Page 5 | **Elbow Connectors**

Page 6 | **Cuff Pressure Monitor**

Page 7 | **Breathing Circuit Tubing**

Page 8 | **Pulse Oximeter Finger Wrap**

Page 9 | **Trach Foam Dressing**

# RESPIRATORY MASKS



**Trach Masks**



**Aerosol Masks**



**Venturi Masks**

**Oxygen Masks**

## Varieties Available:

- **Aerosol Masks** - 360 degree swivel connector for maximum patient comfort
- **Venturi Masks** - Rotating orifice which is able to set up the flow rates and percentage concentration required
- **Oxygen Masks** - Swivel oxygen mask with CO2 sampling line and filter
- **Trach Masks** - 360 degree swivel connector with 6mm connector

## Features and Benefits:

- 100% medical grade material
- Kink-free, crush-resistant tubing
- Clear, soft vinyl for patient comfort and visual patient assessment
- Material: DEHP free, latex-free

### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax  
orders@peppermedical.com | www.PepperMedical.com

## NEBULIZER KITS



### Features and Benefits:

- Disposable Nebulizer Kits
- Mouth pieces
- T-pieces
- Corrugate tubing
- Oxygen tubing
- Jet Nebulizer Sets

#### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax  
orders@peppermedical.com | www.PepperMedical.com

## OXYGEN TUBING



### Features and Benefits:

- 100% medical grade material
- Kink-free, crush-resistant tubing
- Material: DEHP free, latex-free
- Universal connectors on both sides
- Available in various sizes (7 ft, 25 ft, etc.)
- Available with and without cannulas
- Available for adult, pediatric and infant
- Available with CO2 sampling line

#### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax  
orders@peppermedical.com | www.PepperMedical.com

## ELBOW CONNECTORS



### Benefits:

- Disposable Elbow Connectors
- Made from 100% medical grade K resin plastic
- Available with 15mm and 22mm connectors on either or both sides
- Available with and without Luer Locks
- Available with and without oxygen ports
- Swivel options available

### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax  
orders@peppermedical.com | www.PepperMedical.com

# DISPOSABLE CUFF PRESSURE MONITOR



## Description:

Disposable cuff gauge is intended to monitor Tracheostomy Tube and Endotracheal Tube cuff pressures within the safe zone to maintain proper patient care.

## Features and Benefits:

- Easy to use syringe connector for quick cuff pressure adjustment
- Under inflated cuff pressure can cause ventilator leakage and aspiration while over inflated cuff pressures may lead to tracheal stenosis, necrosis, tracheomalacia and long term structural damage to the airway.
- The Pepper Medical Disposable Cuff Pressure Monitor assures clinical evaluation and visibility within a safe and reliable zone. Clinicians no longer need to guess or palpate the pilot balloon.
- The Pepper Medical Disposable Cuff Monitor is a safe, easy and inexpensive way to assure patient comfort and cuff pressure protocols.
- Real Time accurate cuff pressure readings
- Highly visible markings and Safe Zone indicator
- Reduces tracheal mucosal tissue damage
- Inexpensive alternative to expensive gauges
- Single patient use to reduce cross contamination

## Sizes:



- Tracheal- 22~32cmH<sub>2</sub>O; LMA- 42~60cmH<sub>2</sub>O

### ORDERING INFORMATION:


1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

# BREATHING CIRCUIT TUBING



**Corrugated  
Taper Tubing**



**Expandable  
Taper Tubing**

## Benefits:

- Tubing for breathing circuit systems
- Fits standard 15 mm and 22 mm connections
- Available in two configurations:

**Circuit Extension | Expandable, 6"**

**Circuit Extension | Corrugated, 6"**

### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax  
orders@peppermedical.com | www.PepperMedical.com



#654

## PULSE OXIMETER FINGER WRAP



### Benefits:

- Pulse oximeter probe wrap that helps secure various styles of optical transducers and probe sensors
- Prevents the misalignment of optical components while blocking out extra light
- Fits a neonatal foot, pediatric toes, or pediatric and adult fingers
- Packaged 12 wraps per package
- Velcro® tab eliminates need for adhesive tape
- Skin friendly foam laminate made from moisture repellent material
- Over-the-counter intended for non-medical, recreational use
- Adjustable, one size fits all
- Secures Pulse Ox sensor to allow for steady reading
- Measures 5" long x 1" wide and can be cut to fit any length

### ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax  
orders@peppermedical.com | www.PepperMedical.com



## **TRACH FOAM DRESSING**

---

### Benefits:

- Composed of polyurethane foam
  - Suitable for a variety of exudative ostomy wounds, including Tracheostomy
  - Texture is soft, the pore size is uniform which allows for high absorbency
  - Provides the best environment for moist wounds
  - Promotes wound healing
- 

### Features:

- Waterproof
  - High absorbency
  - Water-lock function
  - Air-permeable and breathable
- 

#### **ORDERING INFORMATION:**

1-800-647-0172 | 610-692-8468 Fax  
orders@peppermedical.com | www.PepperMedical.com