

# DISPOSABLE CUFF PRESSURE MONITOR



## Description:

Disposable cuff gauge is intended to monitor Tracheostomy Tube and Endotracheal Tube cuff pressures within the safe zone to maintain proper patient care.

## Features and Benefits:

- Easy to use syringe connector for quick cuff pressure adjustment
- Under inflated cuff pressure can cause ventilator leakage and aspiration while over inflated cuff pressures may lead to tracheal stenosis, necrosis, tracheomalacia and long term structural damage to the airway.
- The Pepper Medical Disposable Cuff Pressure Monitor assures clinical evaluation and visibility within a safe and reliable zone. Clinicians no longer need to guess or palpate the pilot balloon.
- The Pepper Medical Disposable Cuff Monitor is a safe, easy and inexpensive way to assure patient comfort and cuff pressure protocols.
- Real Time accurate cuff pressure readings
- Highly visible markings and Safe Zone indicator
- Reduces tracheal mucosal tissue damage
- Inexpensive alternative to expensive gauges
- Single patient use to reduce cross contamination

## Sizes:

- Tracheal- 22~32cmH<sub>2</sub>O; LMA- 42~60cmH<sub>2</sub>O