



PRODUCT #501W Tracheostomy Tube Neckband, Wide

Features & Benefits:

- Fits up to a 28" neck
- 1.5" wide design prevents roping effect
- Two piece design allows for easy adjustment
- Soft, plush backing on hook tabs provides additional comfort
- Made of DriWick®, a proprietary material that pulls moisture away for the patients skin, reducing skin breakdown and irritation
- LATEX FREE

MOISTURE MANAGEMENT PERFORMANCE DriWick VS Cotton

<u>CATEGORY</u>	<u>PERFORMANCE</u>	
	<u>DriWick</u>	<u>Cotton</u>
WET OUT	INSTANTANEOUS	BEADED –WATER DROP HAD TO BE FORCED INTO FABRIC
DRYING TIME	< 10 MINUTES	> 60 MINUTES
VERTICAL WICKING TEST (INTERNAL NATICK METHOD)	4 INCHES IN 30 MINS	¼ INCH IN 30 MINS