



210 Carter Drive Unit 10, West Chester, PA 19382
Phone 800-647-0172 Fax 610-692-8468 www.peppermedical.com

Item #BV001 – Bacterial Viral Filter



Product Code	#BV001
Description	Bacterial Viral Filter
Weight	G.W.: 5 N.W.: 4
Filtration	BFE 99.9999% VFE 99.9999%
Gas Port	Yes
Connection:	22M-15F / 22F-15M

Manufacturing and General Offices
600 W Strasburg Road • West Chester, PA • 484-653-5473