



600 W Strasburg Road, West Chester, PA 19382

Phone 800-647-0172

Fax 610-692-8468

www.peppermedical.com

HMEF711C



| | |
|---------------------|----------------------------------|
| Product Code | #HMEF711C |
| Description | HME With Filter |
| Humidity Output | 31mg H ₂ O @ VT 500ml |
| Resistance to Flow | 180pa @ 60 LPM |
| Tidal Volume | 150 – 1500 ml |
| Filtration | BFE 99.9999% VFE 99.9999% |
| Dead Space | 45 ml |
| Gas Port | Yes |
| Weight | 30 grams |
| Connection | 22M-15F / 22F-15M |

Manufacturing and General Offices
600 W Strasburg Road • West Chester, PA • 484-653-5473