



**Heat and Moisture Exchanger - Mini
#HME-Mini-P**



Product Code	HME-Mini-P
Description	Mini HME
Humidity Output	32mg H ₂ O @VT 25ml
Resistance @ 5 LPM	90pa
Dead space	3ml
Tidal Volume Range	> 10ml
Weight	3g

The Mini Tracheal Heat & Moisture Exchanger captures heat and moisture and returns it to the patient. It is meant for use by spontaneously breathing patients with a tracheostomy tube.

Product Features

- Low profile and Lightweight.
- Low Volume Tracheal HME
- Connects to a 15F-15M connection