



600 W Strasburg Road, West Chester, PA 19382  
Phone 800-647-0172 Fax 610-692-8468 www.peppermedical.com

### Pepper Medical– #540 series Trach Ties

Pepper Medical Item #	Description
#540 <b>Width:</b> 1" wide <b>Length:</b> 18.5" long including 2.5" of elastic <b>Foam Laminate:</b> 0.25" thick blue foam with moisture repellant backing <b>Tabs:</b> 2 .875" Velcro hook tabs, stitched	2 Piece Trach Tube Holder (Adult)
#541 <b>Width:</b> 1" wide <b>Length:</b> 17.5" long <b>Foam Laminate:</b> 0.1875" thick blue foam with moisture repellant backing <b>Tabs:</b> 2 .875" Velcro hook tabs, welded	2 Piece Trach Tube Holder (Pediatric)
#542 <b>Width:</b> 0.75" wide <b>Length:</b> 10" long <b>Foam Laminate:</b> 0.1875" thick blue foam with moisture repellant backing <b>Tabs:</b> 2 .875" Velcro hook tabs, welded	2 Piece Trach Tube Holder (Neonatal)

