

FOLEY CATHETER HOLDERS

Catalog



PEPPER MEDICAL

1-800-647-0172 | PepperMedical.com

Contents

Page 2 | #**601**


Pages 3 - 4 | #**606 & #606XL**

Page 5 | #**5445**

Page 6 | #**650**

#601

FOLEY SWIVEL LOC® **Foley Catheter** **Stabilization Device**



Features & Benefits:

- Swivel locking mechanism allows for rotation of the Foley Catheter
- Locking system reduces potential of catheter kinking and twisting
- Minimizes dislodgements, traction and positioning
- **LATEX FREE** - No latex reactions which prevents latex hypersensitivity
- Addresses CDC's guidelines for securing Foley Catheters
- Reimbursement Code A4334

DIMENSIONS

Length: 3" / Width: 2"

ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

#606

FOLEY-TIE®

Foley Catheter Legband



Features & Benefits:

- One piece Velcro® blue locking tab “Flex Point”
- Soft stretch leg band material
- Quick and easy application around leg and Foley
- No adhesive material
- **LATEX FREE** - No latex reactions which prevents latex hypersensitivity
- Fits up to 24” (most adults and children)
- Reimbursement Code A4334

DIMENSIONS

Length: 35” / Width: 2”

ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

#606XL

FOLEY-TIE®

Foley Catheter Bariatric Legband



Features & Benefits:

- One piece Velcro® blue locking tab “Flex Point”
- Soft stretch leg band material
- Quick and easy application around leg and Foley
- No adhesive material
- **LATEX FREE** - No latex reactions which prevents latex hypersensitivity
- Fits up to 60”
- Reimbursement Code A4334

DIMENSIONS

Length: 87.5” / Width: 2”

ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

#5445

TUBE-TIE® Patch – Multi-Purpose Adhesive Tube Holder With Quick Release



Applications Include:

- Foley Catheter holder
- Jejunostomy & Gastrostomy tube holders
- Wound drainage tube
- Subclavian line & other ports

Features & Benefits:

- Securely locks catheter in place
- Non-slip Velcro® Locking Tab can be repositioned
- Secures all types of tubing
- Produces quick and effective removals
- **LATEX FREE** - No latex reactions which prevents latex hypersensitivity
- Individually packaged in boxes of 50 each
- Reimbursement Code A4334

DIMENSIONS

Length: 3" / Width: 3"


ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com

#650

Foley Loc® Adhesive Foley Catheter Holder



Features & Benefits:

- Kink Resistant – non-swivel locking system reduces potential of catheter kinking and twisting
- Adjustable – Velcro® EZ-Tabs allow for release, inspection / cleaning and re-attachment
- Security – minimizes dislodgements, traction and pistoning
- Versatility – secures and stabilizes all styles of Foley catheters (silicone or latex, two or three port catheters)
- Universal Design – one style reduces inventory and save costs
- **LATEX FREE** and hypo-allergenic medical grade adhesive
- Addresses CDC's guidelines for securing Foley Catheters
- Individually packaged in boxes of 50 each
- Reimbursement Code A4333

DIMENSIONS

Length: 1.75" / Width: 5.5"

ORDERING INFORMATION:

1-800-647-0172 | 610-692-8468 Fax

orders@peppermedical.com | www.PepperMedical.com