




Speaking Valves *Specifications*

Speaking Valves	
Pepper Medical Item #	Description
<p>#710R – Speaking Valve</p> 	<p>Speaking Valve with oxygen port for spontaneous breathing tracheostomy patients.</p> <p>Intended to allow the patient to speak without the need for finger occlusion. Also provides back pressure for other common body functions. Exchangeable filter to prevent course particles entering the respiratory system.</p>
<p>#710 – Speaking Valve</p> 	<p>Speaking Valve with oxygen port for spontaneous breathing tracheostomy patients. Intended to allow the patient to speak without the need for finger occlusion. Also provides back pressure for other common body functions. Exchangeable filter to prevent course particles entering the respiratory system.</p> <p>This transparent device is recommended for use in the home care environment.</p> <p>*Also available without O2 Port</p>
<p>#740- Replacement Filter</p> 	<p>Exchangeable filter which is compatible 710, 710R. The filter prevents course particles from entering respiratory system.</p>